



Lyophilised Zika Virus detection kit from PCRmax

The Delta Seek detection kits utilise the sensitivity and speed of qPCR to get the most accurate data as quickly as possible. Each kit is specifically designed by our bioinformatics team to ensure the broadest possible detection profile and detection of all clinically relevant strains and subtypes. All test kits are validated in house on multiple qPCR platforms to ensure cross platform functionality.

Delta Seek kits are also customisable. If your target is not included in our extensive range of over 400 kits covering the fields below, we can make custom kits for you.



DeltaSeek qPCR Reagent Kits



Delta seek kits cover the following areas of research:

- Human Pathogens
- Food and Water Testing
- Veterinary and Agriculture
- Biothreat Detection

For a comprehensive and searchable list of all the available Delta Seek kits visit www.PCRmax.com

Features and Benefits:

- Positive copy number standard curve for accurate and absolute quantification.
- Highly specific detection profile with high efficiency priming.
- Rapid detection of all clinically relevant subtypes.
- Sensitive to < 100 copies of target.
- Only two protocols (one DNA and one RNA) covers the full range.
- All contents lyophilised for simple transport and storage.

www.pcrmax.com

PCRmax